



umano medical

Pilot Drive™



A revolutionary powered drive system

Mobility, versatility, patient safety,
and product performance all-in-one.

Reimagined with **ease of use** and **fall prevention** in mind, caregivers can still **benefit from our exceptional low height** and **essential underbed clearance** for lifts and overbed tables.

Infection prevention and control is supported through our renowned seamless smooth surface design – ensuring **cleaning remains an easy task**.

Available on the



Bring out the Pilot in you.

Feel confident and in control when transporting
ook«SNOW» ALL beds



Offering variable speed capabilities and battery autonomy, the **Pilot Drive** allows staff to cover long distances, even up ramps and steep slopes.



The **Pilot Drive's** ergonomic design and powerful driving force helps relieve excessive push and pull stress for staff.

TECHNICAL SPECIFICATIONS

Bed's lowest height	12.7" (32.3 cm)
Underbed clearance for tables and lifts	5" (12.7 cm) at bed's lowest height
Variable travel speed	0 to 3.1 mph (5 km/h)
5th wheel	Centrally located to make tight and controlled turns
Inclination	Maneuverable in slopes of up to 5°
Options	Integrated oxygen bottle holders (compatible with E size bottles)

Umano Medical reserves the right to change specifications without notice.



Umano Medical Inc.
230, boulevard Nilus-Leclerc
L'Islet QC G0R 2C0 CANADA
Toll Free: 1 (844) 409-4030
information@umanomedical.com
umanomedical.com

Bring out the Pilot in you.