



# Integrity 270

## Steam Sterilizers

Integrity 270 is the workhorse of the Integrity line, exceeding all standards for capacity and processing up to sixteen 25-pound trays per cycle. Building on Skytron's history of quality and innovation, the Integrity 270 is part of a family of sterilizers designed to meet the wide-ranging needs of today's healthcare environments. Whether you're outfitting a new facility, renovating an older one, or replacing outdated technology, Integrity sterilizers offer industry-changing efficiency, capacity, ease of use, and use of space.

### **Unprecedented throughput**

Accommodates carriage and loading cart with up to sixteen 25-pound trays in a tight footprint.

### **Save on energy, water, and chemicals**

Integrated water recirculation reduces water usage by up to 67%\*, without affecting cycle time or increasing footprint.

### **Long life and trouble-free operation**

The Integrity Series offers fully jacketed stainless steel chambers, chain-driven electric-powered doors, and stainless steel and copper piping standard.

### **Superb diagnostics and support**

System-wide parts monitoring, on-screen diagnostics, local technicians, and 24/7/365 call center ensure the reliable operation of your equipment.

### **Offers high efficiency evacuation of air and steam**

Optional vacuum pump available.

### **Simplifies replacement of older technology**

Integrity 270's plug-and-play installation makes replacing outdated product easy.



# Integrity 270

## Steam Sterilizers

24 Preset and custom cycles Sterilize all your steam-compatible devices and instrument sets with simple touch-screen control.

### Processing cycles

- Immediate use 1 (3 minute)
- Immediate use 2 (4 minute)
- Express prevacuum (immediate use 3)
- Prevacuum 270
- Prevacuum 275
- Gravity 1
- Gravity 2
- Liquid cycle (optional)
- Preset cycles validated to AAMI ST-8:2008

### Test cycles

- Warm up
- BD
- Cycle
- Leak

### Custom cycles

12 custom cycles modifiable to your unique needs

Freestanding stainless steel cabinet available for single door models.

- Spacious 26½" x 37½" x 53" sterilization chamber
- Processes up to sixteen 25-pound trays (400 pounds per cycle)
- Integrated water recirculation system
- 5.7" LCD full-color touchscreen panel
- Countdown timer and progression bar visible from more than 30'
- 12 standard and 12 fully customizable cycles
- Electric-powered door locking/unlocking mechanism
- Stainless steel and copper piping
- Dual-walled, fully jacketed stainless steel pressure vessel
- Electronic recording of cycle report data
- Integrated ink-to-paper printer for hard copy cycle reports
- Programmed automatic cycle start
- Energy saving daily startup/shutdown timer
- Password protected access to supervisor and service modes
- Compliant to AAMI-ST8:2008 and AAMI-ST79:2010 standards

### Options

- Recessed mounting or freestanding cabinet
- Standard jet pump or optional vacuum pump (VP requires high-voltage electrical)
- Optional liquid cycle and liquid temperature load probe
- Single or double door (DD) that swings right or left, as specified
- Right or left hand service access, as specified

Order Configurations							
Model	Door Type	Vacuum	High Voltage	Mounting	Door Swing	Service Side	Equipment #
Integrity 270	Single	Jet Pump	N/A	Recessed	Right	Right	ASSR-A012-RRSR
						Left	ASSR-A012-RRSL
					Left	Right	ASSR-A012-RLSR
						Left	ASSR-A012-RLSL
	Double			Recessed	Right	ASSR-A012W-RRSR	
					Left	ASSR-A012W-RRSL	
				Left	Right	ASSR-A012W-RLSR	
					Left	ASSR-A012W-RLSL	
Integrity 270VP	Single	Vacuum Pump	208V, 3pH, 60Hz, 15A	Recessed	Right	Right	ASSR-A012P-208R-RSR
						Left	ASSR-A012P-208R-RSL
					Left	Right	ASSR-A012P-208R-LSR
				Left		ASSR-A012P-208R-LSL	
				Recessed	Right	Right	ASSR-A012P-480R-RSR
						Left	ASSR-A012P-480R-RSL
	Left		Right		ASSR-A012P-480R-LSR		
			Left	ASSR-A012P-480R-LSL			
	Double		Recessed	208V, 3pH, 60Hz, 15A	Right	Right	ASSR-A012PW-208R-RSL
						Left	ASSR-A012PW-208R-RSR
				Left	Right	ASSR-A012PW-480R-LSL	
			Left		ASSR-A012PW-480R-LSR		
Recessed		208V, 3pH, 60Hz, 15A	Right	Right	ASSR-A012PW-208R-RSL		
				Left	ASSR-A012PW-208R-RSR		
	Left	Right	ASSR-A012PW-480R-LSL				
Left		ASSR-A012PW-480R-LSR					

\*Water usage is based on comparison of fully loaded sterilizers with similar-size chambers, using a 270°F prevacuum cycle with 4-minute sterilize and 30-minute dry time. Testing showed that the Integrity 270 used 33% less water than the similarly sized STERIS® EVOLUTION™.

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at [www.skytron.com](http://www.skytron.com)



5085 Corporate Exchange Blvd. SE  
 Grand Rapids, MI 49512  
 1.800.SKYTRON (759-8766)  
 email: [info@skytron.com](mailto:info@skytron.com)  
[www.skytron.com](http://www.skytron.com)

