



# GS70 Guardian Series

## G-Series Table Accessories

The G-Series accessories were designed specifically for use on all Skytron G-Series surgical tables. When an accessory is attached to the leg section, the table will automatically recognize the unique dimensions of the accessory and will prevent any collisions with the floor or the table base. Skytron G-Series accessories are only for use on Skytron model G-Series tables and are not compatible with older model tables.

### **Safety for Staff and Patients**

The accessories are made of high-quality materials and are tested and verified for use on Skytron tables.

### **Less Damages, More Savings**

With the G-Series smart technology, the table can alert staff before a potential accident can occur.

### **Easy to Move**

Made with carbon fiber materials, the GS70 accessories are lightweight making it easier to attach and remove accessories from the table.




# GS70 Guardian

## G-Series Table Accessories


### 44" Carbon Fiber Extension - PN 4-031-01

The 44" extension provides unmatched imaging and surgical access all while maintaining an 800 lb. patient capacity.




### Kidney Ureter Bladder - KUB - PN 4-041-01

The Kidney Ureter Bladder board provides additional imaging length and width while the carbon fiber allows for 360° of imaging making oblique c-arm imaging easy.



### Light Weight Leg Section - PN 4-031-02

The lightweight leg section provides an alternative to the standard leg section which is lighter and easier to take on and off. This accessory is most useful for constant lithotomy positioning when the leg section is regularly removed and reinstalled.




### Split Leg - PN 4-030-70

Powered by the table itself the split leg ingrates perfectly with G-Series tables. The split leg works independently or in unison by raising or lowering each leg using table power.




### Pediatric Extension - PN 4-031-03

The pediatric extensions provides direct access to smaller patients. Eliminate offsetting patients or leaning over the table. At 15" wide this extension provides great access from either side for your smallest patients.



### 5" Side Extension - PN 3-060-09

The 5" table side extensions will not slide or lift from the table side rails and lock firmly in place. They will provide a 30" working surface to meet bariatric needs.



Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at [www.skytron.com](http://www.skytron.com)



5085 Corporate Exchange Blvd. SE  
Grand Rapids, MI 49512  
1.800.SKYTRON (759-8766)  
email: [info@skytron.com](mailto:info@skytron.com)  
[www.skytron.com](http://www.skytron.com)



REV1 09-22  
PN L2-010-200