



Processing Sink

Fixed or Height-Adjustable Solution for
Sterile Processing Departments



Sterile Processing Expertise that's Designed with Superior Workflow in Mind

Processing sinks serve as the pivotal hub of the sterile processing department where meticulous cleanliness and care meets efficiency. Skytron's new processing sink complements your disinfection processes and provides a means for the pre-cleaning, inspection, and sorting of surgical instruments.

Configurable to Your Sterile Processing Setting

- Single, dual, or triple basins to fit your department needs best
- Ergonomically friendly design with an electronic height adjustable range of 35"-42"
- Universal Workflow design that complements all SPD working styles
- Optional automated dosing system ensures the correct amount of detergent is used, resulting in long-term cost savings for the facility
- Backsplash-mounted faucets and air and water guns free up work surface around basins
- Optional pegboard for added flexibility to easily configure accessories
- Larger basin option to support robotic instruments and overall increased tool sizes



Fully
Configurable
Menu



Optional Electronic
Height Adjustability



Compatible
with Clean
Start Chemicals

Functional Organizational Accessory Options

Standard Features

- Single, dual, or triple basins
- Universal workflow
- Fixed and height adjustable (37" standard height for fixed height sinks)
- Ergonomically friendly electronic adjustable height range: 35"-42"
- Basin size: 24"W x 10"D x 20"L
- Laser-etched gallon markings
- Marine edge working surface
- Fully welded construction with reinforced tops
- Adjustable feet
- Easy access to connection points with open pipe base

New Standard Features

- Universal workflow
- 12" angled backsplash
- Backsplash mounted faucets
- 30" deep standard sink
- 20" deep standard basin

Optional Features

- Pre-rinse faucet
- Mixing valve faucet
- Critical water faucet
- Air gun kit (with 5 tips)
- Water gun kit (with 8 tips)
- Hand lever drains
- Pegboard and optional accessories (sold separately)*
- Copper drain system**
- Basin size upgrade: 30"W x 10"D x 20"L

New Optional Features

- Pegboard backframe and accessories
- Side splash (bolt-on)
- Automated detergent dosing
- Endoscope cleaning system
- Automated scope flushing system
- Light bar with LED light
- Bin rail with bins

** Pegboard optional accessories include glove dispenser, gallon jug holder, spray bottle holder, overhead shelf, dosing system, and monitor mounts.*

***Copper drain system is required for compliance in CA*



Cleaner Work Environment

Skytron's Processing Sinks come standard with 20" deep basins. The basins are laser-etched in 5 gallon increments with no crevices that could house harmful bacteria.



Alleviate Physical Stress

The height adjustable option ranging from 35"-42" helps alleviate any physical stress.



Clean Start Chemicals

Super concentrated for cost effectiveness. Extends instrument life with anti-corrosive agents.

Skytron's Solutions



Sterile Processing

Processing Sinks and Stainless Offerings
Washers, Ultrasonics, and Cart Washers
Prep & Pack Tables & Sterilizers



Clinical

Surgical Tables
Booms & Lights
Full Offering of Stainless Products



Architectural & Environmental

Modular Walls
Configurable Window and Door Offerings
Built-in Cabinetry



Clinical Business Intelligence

Video Integration
Asset Tracking & Supply Chain Solutions
Surgical Displays

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service. Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
Grand Rapids, MI 49512
1.800.SKYTRON (759-8766)
email: info@skytron.com
www.skytron.com



REV1 06-24
PN L2-010-256